




### Our mission

We are committed to making a positive difference in the lives of dental practitioners and patients worldwide by bringing products to the market that make dentistry easier and safer. We strive to provide products that enable dentists to fulfill a social responsibility to provide better, healthier patient care. We seek to provide a greater sense of satisfaction to each of our global customers by offering truly innovative, effective, high-value products and solutions supported by excellent quality and customer service.

# NEOSHINE®



### CONTACT US

Sales@microcopyintl.com  
Microcopyintl.com



Microcopy

3120 Moon Station Rd. NW  
Kennesaw, GA 30144, U.S.A.

# NEOSHINE®

## Complete Polishing System

Trouble polishing on the anterior? Thin single-patient-use, NEOSHINE Discs provide flexibility to hug the curves of anterior surfaces. Problem solved! NEOSHINE Discs are available in four grits (very coarse, coarse, medium, and fine) to remove irregularities with ease. Expect a smooth hi-gloss polish on composite and porcelain. With one side being abrasive, NEOSHINE Discs are also able to polish interproximally.

- Flexible disc enables you to thoroughly polish all tooth anatomy.
- Polishes to a smooth mirror-like finish.
- Single sided discs allow you to polish interproximally.

## Discs


Available in two sizes, 10mm and 14mm. One size per pack of 100 discs (25 pcs. per grit) with one reusable stainless-steel RA mandrel.

## Instructions for Use

- **Composite/Porcelain:** Very Coarse, Coarse, Medium & Fine.
- **Recommended Speed:** 10,000 rpm **Max speed:** 20,000 rpm.\*
- Polishing paste is not required.
- Keep the disc dry and cool tooth with air at least every 10 to 15 seconds to avoid excess heat build-up.

*\*Avoid long periods of applied pressure as heat generation can cause damage to the tooth. Improper use of these instruments may lead to poor results and increased patient risk.*

## COMPOSITE / PORCELAIN - DISC



Very Coarse Coarse Medium Fine

### Order No.

(10mm)  
NS-CP DISC10

(14mm)  
NS-CP DISC14

One size per pack of 100 pcs.  
(25 pcs. per grit) + 1 mandrel

## DISC REFILL



Very Coarse  
NS-CPVC DISC 10  
NS-CPVC DISC 14



Coarse  
NS-CPC DISC 10  
NS-CPC DISC 14



Medium  
NS-CPM DISC 10  
NS-CPM DISC 14



Fine  
NS-CPF DISC 10  
NS-CPF DISC 14

### Order No.

See image above

One size and grit per pack  
of 100 pcs.


## RA MANDREL REFILL



3 per pack

"The speed of polishing  
and quality of the resulting surface  
finish are second to none!"

Joshua Howard DMD  
Dover, NH



⊗ Do no re-use ⓘ nspifu.com