


NEODIAMOND® Kits



Crown Preparation




Cavity Preparation




Endo



Composite & Porcelain Finishing




Laminate Veneer




Pediatric



Zirconia Finishing



NeoSpiral



NEODIAMOND® Application Chart

Round 001 801 801L Crown preparation	Round with Collar 002 802 Crown preparation	Inverted Cone short head 010 805 Removal of old fillings Cavity preparation	Inverted Cone with collar 019 806 Removal of old fillings Cavity preparation	Diablo 032 813 Removal of old fillings Cavity preparation	Barrel 038 039 811 Crown preparation Occlusal contouring Contouring+Finishing
--	--	---	--	---	--



Wheel 041 067 068 815 909 Crown preparation Occlusal contouring	Wheel with shoulder 055 824 Cavity preparation	Flat End Cylinder 108 109 110 111 112 835 836 837 Removal of old fillings Cavity preparation Crown preparation Inlay/Onlay preparation	Beveled Cylinder 129 130 131 884 885 886 Crown preparation Removal of old fillings	Round End Cylinder short head 138 139 838 Removal of old fillings Cavity preparation Crown preparation Ceramic Restorations Inlay/Onlay preparation
--	---	--	--	--

Round End Cylinder 140 141 880 881 Removal of old fillings Cavity preparation Crown preparation Ceramic Restorations Inlay/Onlay preparation	End Cutting Cylinder 150 839 Cavity preparation Crown preparation	Round Edge Cylinder 164 165 166 167 852 858 859 Removal of old fillings Cavity preparation Crown preparation Ceramic Restorations Inlay/Onlay preparation	Pointed Cone 164 165 166 167 852 858 859 Crown preparation Surface Finishing	Flat End Taper short head 169 170 845 Cavity preparation Inlay/Onlay preparation
--	---	--	--	--

Flat End Taper 171 172 173 174 846 Removal of old fillings Cavity preparation Crown preparation Inlay/Onlay preparation	Round End Taper 196 197 198 199 200 849 855 856 850 Removal of old fillings Crown preparation Laminate Preparation	Inverted Cone 225 226 233 234 807 830 830L Cavity preparation Occlusal contouring	Pear 237 238 830R 830RL Removal of old fillings Cavity preparation Inlay/Onlay preparation	Flame 248 249 250 860 861 862 Crown preparation Inlay/Onlay preparation
---	--	---	--	---

Bud 257 368 Occlusal contouring Crown preparation Contouring & Finishing	Egg 277 379 Crown preparation Contouring & Finishing Occlusal contouring	Torpedo 288 289 290 877 878 879 Removal of old fillings Crown preparation	Gingival Curettage 297 297 299 877K 878K 879 Cavity preparation Crown preparation	Lenticular 304 825 Crown preparation
---	---	---	---	---

Interdental 465 392 Crown preparation Surface Finishing	Guide Pin 508C 508T Crown preparation	Gross Reduction Cylinder 511 512 513 514 Taper 516 517 519 521 Crown preparation	Flame Needle short head 536 539 540 883 883 889 Crown preparation Surface Finishing	Round Edge Taper 544 545 546 552 834 845 KR 846 KR 847 KR Cavity preparation Inlay/Onlay preparation	Depth Marker 552 834 Crown preparation Laminate Preparation
---	---	---	---	--	---